

EMC-Pad for UV/VIS Spectrophotometer

- Intel® Atom™ Processor x5-Z8350 (2M Cache, up to 1.92 GHz)
- Windows® 10
- Touch Display 10.8" (27.4 cm)
- 64 GB (EMMC) Flash-memory
- 4 GB RAM
- Installed PROFESSIONAL Software incl. dongle
- Stand for EMC-Pad
- USB hub 4-port



EMC-Pad without pictured EMC-Spectrophotometer!

Complete set EMC-Pad with PROFESSIONAL Software and accessories.

The PROFESSIONAL Software is optionally available for Spectrophotometers EMC-11-UV as well as EMC-16/18 series.

Technical changes reserved